




Orthodontic Products




SINGAPORE SMART TECH PTE LTD

U-Click[®] Self Ligating Bracket

The unique sliding clip is made of high quality NiTi, with remarkable anti-fatigue and anti-deformation capabilities, it provides excellent sliding mechanics and delivers consistent tooth moving force, ensures fast tooth movement while maintaining maximum control.

Available in monoblock or with mesh base.



U Click Self Ligating Bracket-Monoblock

		slot	w/h on 3	w/h on 345
Roth	1 case	.022	UE301-020	UE301-030
		.018	UE301-120	UE301-130
MBT	1 case	.022	UE302-020	UE302-030
		.018	UE302-120	UE302-130


U Click Self Ligating Bracket-Mesh Base

		slot	w/h on 3	w/h on 345
Roth	1 case	.022	UM301-020	UM301-030
		.018	UM301-120	UM301-130
MBT	1 case	.022	UM302-020	UM302-030
		.018	UM302-120	UM302-130

Sapphire Bracket

The bracket is created from a single crystal of pure grown sapphire – one of the hardest materials in nature, second only to diamond. The crystal is honed and heat polished into a remarkably strong and beautiful bracket which is nearly invisible on the teeth. The crystal-clear sapphire provides excellent aesthetics. However, the aesthetics can create challenges for a clinician to accurately position brackets. The Surface Paint technology is designed to tackle this very problem. Tie-wings of the brackets are coated with a water-soluble dye creating a crosshair effect that enables easy and precise placement which further enhances the treatment outcome and avoids bracket repositioning, thus saving time.

- Crystal-clear and nearly invisible.
- Excellent ductility.
- Unique Face Painting(FDA edible colorant, water-soluble), colorful brackets easy to be distinguished and positioned.



Sapphire Bracket


		slot	w/h on 3	w/h on 345
Roth	1 case	.022	SP201-020	SP201-030
		.018	SP201-120	SP201-130
MBT	1 case	.022	SP202-020	SP202-030
		.018	SP202-120	SP202-130

Metal Bracket



Metal Bracket-Mesh Base

		slot	no hook	w/h on 3	w/h on 3,4,5
Roth	1 case	.022	BR001-010	BR001-020	BR001-030
		.018	BR001-110	BR001-120	BR001-130
MBT	1 case	.022	BR002-010	BR002-020	BR002-030
		.018	BR002-110	BR002-120	BR002-130



Metal Bracket-Monoblock


		slot	no hook	w/h on 3	w/h on 3,4,5
Roth	1 case	.022	BR003-010	BR003-020	BR003-030
		.018	BR003-110	BR003-120	BR003-130
MBT	1 case	.022	BR004-010	BR004-020	BR004-030
		.018	BR004-110	BR004-120	BR004-130

Ceramic Bracket



Ceramic Bracket


		slot	w/h on 3	w/h on 3,4,5
Roth	1 case	.022	BR201-020	BR201-030
		.018	BR201-120	BR201-130
MBT	1 case	.022	BR202-020	BR202-030
		.018	BR202-120	BR202-130



Ceramic Bracket With Metal Slot

		slot	w/h on 3	w/h on 3,4,5
Roth	1 case	.022	BR301-020	BR301-030
		.018	BR301-120	BR301-130
MBT	1 case	.022	BR302-020	BR302-030
		.018	BR302-120	BR302-130

Buccal Tubes for First Molar



Weldable Single Non-convertible

	slot	UR	UL	LL	LR
Roth	.022	BT001-016	BT001-026	BT001-036	BT001-046
	.018	BT001-116	BT001-126	BT001-136	BT001-146
MBT	.022	BT002-016	BT002-026	BT002-036	BT002-046
	.018	BT002-116	BT002-126	BT002-136	BT002-146
Edgewise	.022	BT003-016	BT003-026	BT003-016	BT003-026
	.018	BT003-116	BT003-126	BT003-116	BT003-126



Weldable Double Rectangular Non-convertible

	slot	UR	UL	LL	LR
Roth	.022	BT005-016	BT005-026	BT005-036	BT005-046
	.018	BT005-116	BT005-126	BT005-136	BT005-146
MBT	.022	BT006-016	BT006-026	BT006-036	BT006-046
	.018	BT006-116	BT006-126	BT006-136	BT006-146
Edgewise	.022	BT007-016	BT007-026	BT007-016	BT007-026
	.018	BT007-116	BT007-126	BT007-116	BT007-126



Weldable Double Combination Non-convertible

	slot	UR	UL	LL	LR
Roth	.022	BT008-016	BT008-026	BT008-036	BT008-046
	.018	BT008-116	BT008-126	BT008-136	BT008-146
MBT	.022	BT009-016	BT009-026	BT009-036	BT009-046
	.018	BT009-116	BT009-126	BT009-136	BT009-146
Edgewise	.022	BT010-016	BT010-026	BT010-036	BT010-046
	.018	BT010-116	BT010-126	BT010-136	BT010-146

Buccal Tubes for First Molar

Weldable Single Convertible



	slot	UR	UL	LL	LR
Roth	.022	BT051-016	BT051-026	BT051-036	BT051-046
	.018	BT051-116	BT051-126	BT051-136	BT051-146
MBT	.022	BT052-016	BT052-026	BT052-036	BT052-046
	.018	BT052-116	BT052-126	BT052-136	BT052-146
Edgewise	.022	BT053-016	BT053-026	BT053-016	BT053-026
	.018	BT053-116	BT053-126	BT053-116	BT053-126

Weldable Double Rectangular Convertible




	slot	UR	UL	LL	LR
Roth	.022	BT054-016	BT054-026	BT054-036	BT054-046
	.018	BT054-116	BT054-126	BT054-136	BT054-146
MBT	.022	BT055-016	BT055-026	BT055-036	BT055-046
	.018	BT055-116	BT055-126	BT055-136	BT055-146
Edgewise	.022	BT056-016	BT056-026	BT056-016	BT056-026
	.018	BT056-116	BT056-126	BT056-116	BT056-126

Weldable Double Combination Convertible



	slot	UR	UL	LL	LR
Roth	.022	BT057-016	BT057-026	BT057-036	BT057-046
	.018	BT057-116	BT057-126	BT057-136	BT057-146
MBT	.022	BT058-016	BT058-026	BT058-036	BT058-046
	.018	BT058-116	BT058-126	BT058-136	BT058-146
Edgewise	.022	BT059-016	BT059-026	BT059-036	BT059-046
	.018	BT059-116	BT059-126	BT059-136	BT059-146

Buccal Tubes for First Molar



Bondable Single Non-convertible

	slot	UR	UL	LL	LR
Roth	.022	BT101-016	BT101-026	BT101-036	BT101-046
	.018	BT101-116	BT101-126	BT101-136	BT101-146
MBT	.022	BT102-016	BT102-026	BT102-036	BT102-046
	.018	BT102-116	BT102-126	BT102-136	BT102-146
Edgewise	.022	BT103-016	BT103-026	BT103-016	BT103-026
	.018	BT103-116	BT103-126	BT103-116	BT103-126



Bondable Double Rectangular Non-convertible

	slot	UR	UL	LL	LR
Roth	.022	BT105-016	BT105-026	BT105-036	BT105-046
	.018	BT105-116	BT105-126	BT105-136	BT105-146
MBT	.022	BT106-016	BT106-026	BT106-036	BT106-046
	.018	BT106-116	BT106-126	BT106-136	BT106-146
Edgewise	.022	BT107-016	BT107-026	BT107-016	BT107-026
	.018	BT107-116	BT107-126	BT107-116	BT107-126



Bondable Double Combination Non-convertible

	slot	UR	UL	LL	LR
Roth	.022	BT108-016	BT108-026	BT108-036	BT108-046
	.018	BT108-116	BT108-126	BT108-136	BT108-146
MBT	.022	BT109-016	BT109-026	BT109-036	BT109-046
	.018	BT109-116	BT109-126	BT109-136	BT109-146
Edgewise	.022	BT110-016	BT110-026	BT110-036	BT110-046
	.018	BT110-116	BT110-126	BT110-136	BT110-146

Buccal Tubes for First Molar



Bondable Single Convertible

	slot	UR	UL	LL	LR
Roth	.022	BT151-016	BT151-026	BT151-036	BT151-046
	.018	BT151-116	BT151-126	BT151-136	BT151-146
MBT	.022	BT152-016	BT152-026	BT152-036	BT152-046
	.018	BT152-116	BT152-126	BT152-136	BT152-146
Edgewise	.022	BT153-016	BT153-026	BT153-016	BT153-026
	.018	BT153-116	BT153-126	BT153-116	BT153-126



Bondable Double Rectangular Convertible

	slot	UR	UL	LL	LR
Roth	.022	BT154-016	BT154-026	BT154-036	BT154-046
	.018	BT154-116	BT154-126	BT154-136	BT154-146
MBT	.022	BT155-016	BT155-026	BT155-036	BT155-046
	.018	BT155-116	BT155-126	BT155-136	BT155-146
Edgewise	.022	BT156-016	BT156-026	BT156-016	BT156-026
	.018	BT156-116	BT156-126	BT156-116	BT156-126



Bondable Double Combination Convertible

	slot	UR	UL	LL	LR
Roth	.022	BT157-016	BT157-026	BT157-036	BT157-046
	.018	BT157-116	BT157-126	BT157-136	BT157-146
MBT	.022	BT158-016	BT158-026	BT158-036	BT158-046
	.018	BT158-116	BT158-126	BT158-136	BT158-146
Edgewise	.022	BT159-016	BT159-026	BT159-036	BT159-046
	.018	BT159-116	BT159-126	BT159-136	BT159-146

Buccal Tubes for Second Molar



Weldable Single Non-Convertible

	slot	UR	UL	LL	LR
Roth	.022	BT001-017	BT001-027	BT001-037	BT001-047
	.018	BT001-117	BT001-127	BT001-137	BT001-147
MBT	.022	BT002-017	BT002-027	BT002-037	BT002-047
	.018	BT002-117	BT002-127	BT002-137	BT002-147
Edgewise	.022	BT003-017	BT003-027	BT003-017	BT003-027
	.018	BT003-117	BT003-127	BT003-117	BT003-127



Weldable Single Convertible

	slot	UR	UL	LL	LR
Roth	.022	BT051-017	BT051-027	BT051-037	BT051-047
	.018	BT051-117	BT051-127	BT051-137	BT051-147
MBT	.022	BT052-017	BT052-027	BT052-037	BT052-047
	.018	BT052-117	BT052-127	BT052-137	BT052-147
Edgewise	.022	BT053-017	BT053-027	BT053-017	BT053-027
	.018	BT053-117	BT053-127	BT053-117	BT053-127



Bondable Single Non-Convertible

	slot	UR	UL	LL	LR
Roth	.022	BT101-017	BT101-027	BT101-037	BT101-047
	.018	BT101-117	BT101-127	BT101-137	BT101-147
MBT	.022	BT102-017	BT102-027	BT102-037	BT102-047
	.018	BT102-117	BT102-127	BT102-137	BT102-147
Edgewise	.022	BT103-017	BT103-027	BT103-017	BT103-027
	.018	BT103-117	BT103-127	BT103-117	BT103-127



Bondable Single Convertible

	slot	UR	UL	LL	LR
Roth	.022	BT151-017	BT151-027	BT151-037	BT151-047
	.018	BT151-117	BT151-127	BT151-137	BT151-147
MBT	.022	BT152-017	BT152-027	BT152-037	BT152-047
	.018	BT152-117	BT152-127	BT152-137	BT152-147
Edgewise	.022	BT153-017	BT153-027	BT153-017	BT153-027
	.018	BT153-117	BT153-127	BT153-117	BT153-127

Pictures



First Molar Bands

Size	Circumference(mm)
28	28.0
28+	28.5
29	29.0
29+	29.5
30	30.0
30+	30.5
31	31.0
31+	31.5
32	32.0
32+	32.5
33	33.0
33+	33.5
34	34.0
34+	34.5
35	35.0
35+	35.5
36	36.0
36+	36.5
37	37.0
37+	37.5
38	38.0
38+	38.5
39	39.0
39+	39.5
40	40.0
40+	40.5
41	41.0
41+	41.5
42	42.0
42+	42.5
43	43.0
43+	43.5
44	44.0

Bands : Sandblasted

 <p>Lingual Buttons 10pcs/bag Bondable Round Base IA04-051</p>	 <p>Lingual Cleats 10pcs/bag Weldable IA10-001</p>								
 <p>Short Crimpable Hooks 10pcs/bag IA05-001</p>	 <p>Lingual Cleats 10pcs/bag Bondable Round Base IA10-051</p>								
 <p>Long Crimpable Hooks 10pcs/bag IA05-101</p>	<p>Ligature Wires</p> <table data-bbox="1084 831 1365 995"> <tbody> <tr> <td>Diameter</td> <td>100g</td> </tr> <tr> <td>0.20m m</td> <td>IA01-120</td> </tr> <tr> <td>0.25m m</td> <td>IA01-125</td> </tr> <tr> <td>0.30m m</td> <td>IA01-130</td> </tr> </tbody> </table>	Diameter	100g	0.20m m	IA01-120	0.25m m	IA01-125	0.30m m	IA01-130
Diameter	100g								
0.20m m	IA01-120								
0.25m m	IA01-125								
0.30m m	IA01-130								

SINGAPORE SMART TECH PTE LTD

Address: 21 Bukit Batok Crescent #15-80 Wcega Tower, Singapore 658065

Tel: 65-66847928 Mobile: 65-92995891 Fax: 65-66847921

Email: sales@singaporesent.com.sg